



Enviral Tech
Zone Surface Check Kit
Dual pathogen: COVID and influenza

Viral Testing Report

COVID-19 & influenza A/B



County Name

Number of COVID-19 cases week of 11/09 **XX**

Number of Influenza cases week of 11/09 **XX**

Building Information		NAME OF COMPANY HERE			
Zone tested	Surfaces swabbed	Collection Date	Test Date	COVID-19 Results	Influenza Results
Zone 1 General Common Area	Elevator buttons	11/12/20	11/13/20	Negative	Negative
	Door handle				
	Table 1				
	Table 2				
	Bathroom door				
	Toilet flush				
	Chair armrest				
	Desk				
Zone 2 Staff Room	Nursing cart 1	11/12/20	11/13/20	Positive	Negative
	Nursing cart 2				
	Employee clock				
	Bathroom door				
	Bathroom faucet				
	Toilet flush				
	Computer keyboard				
	Computer chair				

Legal disclaimer: This report does not constitute legal advice or clinical representation of the presence or absence of the virus COVID-19 or influenza. This report cannot be used to rule out the presence of virus in a facility. It provides only the results from the zones tested. Positive test results indicate the presence of viral particles and may indicate that the zone was contaminated by contact or proximity with an infected person. Enviral Tech encourages individuals and organizations to take universal precautions based on the most current advice from public health officials. For additional information on COVID-19, please see <https://www.cdc.gov/coronavirus/2019-nCoV/index.html>.

Explanation of Results

Negative No Virus (SARS-CoV-2 or influenza A/B) was detected

Positive PCR Testing shows presence of virus in customer sample
Low Positive – between 100 - 200 viral particles
Medium Positive - between 200 - 1000 viral particles
High Positive - >1000 viral particles

Inconclusive The presence of virus (SARS-CoV-2 or influenza A/B) could not be ruled out because of one of these factors:

- Sample was not taken
- RT-PCR reaction was adversely affected by chemicals from the environmental surface that interfere with detection
- A very low signal was seen, but it was below the detection limit of the assay

Note: A positive result is given only if all of the following is true:

1. The positive control tests positive,
2. The negative control tests negative,
3. There is a confirmed sampling using our internal control, and
4. The test of the sample shows the presence of virus

- If either control sample fails, all tests are invalidated and Enviral Tech will rerun the samples.
- If the internal control failed, Enviral Tech will rerun the sample.
- If the internal control fails a second time, the conclusion is that no sample was taken, or that the sample includes environmental contaminants that interfere with the test. We would recommend ordering a new sample kit to retest.